

S/5™ Aespire

Essential performance
Compact design



Shown with Cardiocap™/5 monitor and Tec 7 Vaporizers

Features

- Enhanced monitor integration capabilities with our S/5™ Anesthesia Monitor and Cardiocap™/5 monitor
- Lightweight and compact for easy maneuverability
- Optional integrated auxiliary O₂ flowmeter and suction control

Advanced Breathing System (ABS)

- Bag/vent switch turns the ventilator on/off
- Minimal number of parts and tube connections greatly reduces the potential for leaks and misconnects
- Ease of disassembly (no tools)
- Fully autoclavable and latex-free

Ventilation

- Volume and Pressure modes with electronic PEEP
- Exhaled volume, airway pressure and inspired oxygen monitoring capabilities
- Direct access to ventilator parameter settings
- Smart alarms direct user to specific problems and affected parameters
- Pressure bar graph for visual reference on a breath-by-breath basis (optional pressure waveform available)

Improved low flow/reduced life cycle costs

- One scheduled maintenance check per year
- Fresh gas flow compensation—automatically (available with volume ventilation option)
- Minimum O₂ flow of 50 mL (available with dual flow tube option)
- Optional dual air flow tubes for resolution of low flows

Physical specifications

Dimensions

Height:	134.5 cm/52.9 in
Width:	72 cm/28.3 in
Depth:	73 cm/28.7 in
Weight:	approximately 108 kg/238 lbs

Top shelf

Weight limit:	34 kg/75 lbs
Width:	66 cm/26 in
Depth:	40 cm/15.75 in

Work surface

Height:	81.7 cm/32.2 in
Size:	2640 cm ² /409 in ²

DIN rail

Side of machine:	34.5 cm/13.6 in
------------------	-----------------

Drawers (internal dimensions)



Height:	17.5 cm/6.9 in
Width:	33 cm/13 in
Depth:	26.5 cm/10.4 in

Absorber bag arm (optional)


Arm length:	30.5 cm/12 in
Bag arm height (adjustable):	87 cm/34.3 in 104 cm/40.9 in

Casters

Diameter:	10.1 cm/4 in
Brakes:	Individual locking front casters



Ventilator operating specifications



Optional pressure waveform shown

Modes of ventilation

Volume Control mode

With tidal volume compensation (optional)

Pressure mode (optional)

Ventilation parameters

Tidal volume range: 45 to 1500 mL (Volume Control mode)

Incremental

settings: 45 to 100 mL (increments of 5 mL)
100 to 300 mL (increments of 10 mL)
300 to 1000 mL (increments of 25 mL)
1000 to 1500 mL (increments of 50 mL)

Pressure

(P_{inspired}) range: 5 to 50 cm H₂O (increments of 1 cm H₂O)
5 to 1500 mL volume delivery

Rate: 4 to 65 breaths per minute
(increments of 1 breath per minute)

Inspiratory/
expiratory ratio: 2:1 to 1:6 (increments of 0.5)

Inspiratory
pause adjust: 5% to 60% of inspiratory time (increments of 5%)

Positive End Expiratory Pressure (PEEP)

Type: Integrated, electronically controlled

Range: OFF, 4 to 30 cm H₂O (increments of 1 cm H₂O)

Ventilator monitored values

Tidal volume: 5 to 1500 mL, 1 mL resolution

Minute volume: 0 to 99.9 L/min, 0.1 L/min resolution

Breathing rate: 0 to 65 breaths per minute,
1 breath per minute resolution

Oxygen
percentage: 5% to 110%, 1% resolution

Airway pressure: -9 to 99 cm H₂O, 1 cm H₂O resolution

Alarm settings

Tidal volume (V_{TE}): Low: OFF, 5 to 1500 mL
High: 20 to 1600 mL, OFF

Minute volume (V_{E}): Low: OFF, 0.1 to 10 L/min
High: 0.5 to 30 L/min, OFF

Inspired oxygen
(FiO_2): Low: 18 to 100%
High: 21 to 100%, OFF

Apnea alarm: Mechanical ventilation ON:
< 5 mL breath measured in 30 seconds
Mechanical ventilation OFF:
< 25 mL breath measured in 30 seconds

Low airway
pressure: Change of < 4 cm H₂O above PEEP

Pressure
(P_{limit}) range: 12 to 99 cm H₂O (increments of 1 cm H₂O)

Sustained airway
pressure: 6 to 30 cm H₂O + PEEP (adjusted based on
ventilator settings)

Subatmospheric
pressure: Paw < -10 cm H₂O

Alarm silence
countdown timer: 120 to 0 seconds

Ventilator accuracy

Delivery/monitoring accuracy

Volume delivery:	> 200 mL = better than $\pm 10\%$ Set TV 75 to 200 mL = better than ± 20 mL < 75 mL = better than ± 15 mL
Pressure (P_{Inspired}) delivery repeatability:	± 2 cm H ₂ O
PEEP delivery repeatability:	± 2 cm H ₂ O
Volume monitoring:	> 200 mL = better than $\pm 10\%$ 75 to 200 mL = better than ± 20 mL < 75 mL = better than ± 15 mL
Pressure monitoring:	Better than ± 2 cm H ₂ O and $\pm 5\%$ of reading (whichever is greater)

Ventilator components

Flow transducer

Type:	Variable orifice flow sensor
Dimensions:	22 mm OD and 15 mm ID
Location:	Inspiratory outlet and expiratory inlet

(Optional autoclavable sensor available)

Oxygen sensor

Type:	Galvanic fuel cell
-------	--------------------

Ventilator pneumatics

Pressure range at inlet:	240 kPa to 700 kPa/35 psig to 100 psig
Peak gas flow:	70 L/min + fresh gas flow
Flow range:	2 to 70 L/min
Flow compensation range:	200 mL/min to 15 L/min

Ventilator screen

Display size:	120 mm x 92 mm
Display density:	1/4 VGA standard

Battery back-up

Backup power:	Demonstrated battery time under typical operating conditions is 90 +minutes when fully charged. Battery time under extreme conditions is 30 minutes.
Battery type:	Internal rechargeable sealed lead acid

Communication port

Serial interface:	Isolated RS-232C compatible port
-------------------	----------------------------------


Anesthetic agent delivery

Delivery

Vaporizers:	Tec 5, Tec 6 Plus, Tec 7
Number of positions:	2
Mounting:	Tool-free installation Selectatec® manifold interlocks and isolates vaporizers



Tec 6 Plus



Tec 7

Electrical specifications

Current leakage

100/120 V: < 300µA

220/240 V: < 500µA

Power

Power input: 100-120 Vac, 50/60 Hz
220-240 Vac, 50/60 Hz

Power cord: Length: 5 m/16.4 ft
Rating: 10A @ 220 Vac or 15A @ 120 Vac

Inlet/outlet modules

	220-240 V	120 V	100 V
System circuit breakers:	8A	15A	15A
Outlets (optional):	4 outlets on back, 3-1A, 1-2A individual breakers, optional isolation transformer	4 outlets on back, 3-2A, 1-3A individual breakers, optional isolation transformer	3 outlets on back, 2-2A, 1-4A individual breakers, optional isolation transformer

Pneumatic specifications

Auxiliary common gas outlet

Connector: ISO 22 mm OD and 15 mm ID

Gas supply

Pipeline input range: 240 kPa to 600 kPa/35 psig to 88 psig

Pipeline connections: DISS-male, DISS-female, DIN 13252, AS4059, F90-116, PrEN737-6, or NIST (ISO 5359)
All fittings available for O₂, N₂O, and Air, and contain pipeline filter and check valve.

Cylinder input: Pin indexed in accordance with CGA-V-1 or DIN (nut and gland); contains input filter and check valve

Note: Maximum 3 cylinders; two inboard mounted, one outboard mounted.

Primary regulator diaphragm minimum burst pressure: 2758 kPa/400 psig

Primary regulator nominal output: ≤ 338 kPa/49 psig
Pin indexed cylinder connections
≤ 407 kPa/59 psig
DIN cylinder connections

O₂ controls

Method: Proportionate decrease of N₂O with reduction in O₂ pressure

Supply failure alarm: Range: 193 kPa to 221 kPa/28 psig to 32 psig
Sounds at maximum volume every 10 seconds

O₂ flush: Range: 35 to 50 L/min

Flowmeters

O₂ ranges: One tube option: 0.20 to 12 L/min;
Minimum O₂ flow: 200 mL/min
Two tubes: 0.05 to 0.95 L/min and 1.0 to 15.0 L/min;
Minimum O₂ flow: 50 mL/min ±25 mL

N₂O ranges: Two tubes: 0 to 0.95 L/min and 1.0 to 10.0 L/min

Air range: One tube option: 1 to 15 L/min

Two tube option: 0 to 0.95 and 1 to 15 L/min (low flow tube optional)

Calibration:	Percent of full scale flow	Accuracy (% of flowrate)
	100	±2.5%
	90	±2.5%
	80	±2.6%
	70	±2.7%
	60	±2.9%
	50	±3.1%
	40	±3.4%
	30	±4.0%
	20	±5.0%
	10	±8.1%

Calibration conditions:* 20°C/68°F
101.3 kPa/760 mmHg

* Different breathing circuit pressures, barometric pressures or temperatures change flowtube accuracy.

Hypoxic guard system

Type: Mechanical Link-25™

Range: Provides a nominal minimum 25% concentration of oxygen in O₂/N₂O mixture

Materials

All materials in contact with patient breathing gases are free of natural rubber latex.

Environmental specifications**System operation**

Temperature:	10° to 40°C/50° to 104°F
Humidity:	15 to 95% relative humidity (non-condensing) per IEC 68-2-3
Altitude:	-440 to 3565 m/500 to 800 mmHg

System storage

Temperature:	-15° to 50°C/-5° to 122°F
Humidity:	10 to 95% relative humidity (including condensing) per IEC 68-2-3
Altitude:	-440 to 5860 m/375 to 800 mmHg
Oxygen cell storage:	-15° to 50°C/5° to 122°F 10 to 95% relative humidity 500 to 800 mmHg

Electromagnetic compatibility

Immunity:	Complies with all requirements of EN 60601-1-2
Emissions:	CISPR 11 group 1 class B
Approvals:	UL 2601-1, CSA C22.2 #601.1 EN/IEC 60601-1 CE 0197

Breathing circuit specifications**Operational modes**

Breathing circuit is circle mode only

Carbon dioxide absorbent canister

Absorbent capacity:	800 g
Integrated expiratory limb water reservoir	

Ports and connectors

Exhalation:	22 mm OD ISO 15 mm ID taper
Inhalation:	22 mm OD ISO 15 mm ID taper
Bag port:	22 mm OD

Pressure gauge

Scale range:	0 to 10 kPa/-20 to 100 cm H ₂ O
--------------	--------------------------------------------

Bag-to-Ventilator switch

Type:	Bi-stable
Control:	Controls ventilator and direction of breathing gas within the circuit

Integrated Adjustable Pressure Limiting (APL) valve

Range:	0.8 to 70 cm H ₂ O
Tactile knob indication at:	30 cm H ₂ O and above
Adjustment range of rotation:	0.8 to 30 cm H ₂ O (0 to 230°) 30 to 70 cm H ₂ O (230 to 330°)

Materials

All materials in contact with exhaled patient gases are autoclavable, except disposable flow sensors and O₂ cell. (Autoclavable flow sensors optional)

All materials in contact with patient gas are free of natural rubber latex.

Breathing circuit specifications, continued

Breathing circuit parameters

Compliance:	Bag mode:	1.82 mL/cm H ₂ O	
	Mechanical mode:	Automatically compensates for compression losses within the absorber and bellows assembly	
Circuit volume:	2.7 L Vent Mode		
	1.2 L Bag Mode		
Expiratory resistance:		P _{exp} Bag Mode	P _{exp} Vent Mode
	Flow rate	Pressure drop	Pressure drop
	10 L/min	0.78 cm H ₂ O	0.77 cm H ₂ O
	30 L/min	1.59 cm H ₂ O	1.71 cm H ₂ O
	60 L/min	3.48 cm H ₂ O	3.88 cm H ₂ O

Note: With patient circuit and wye piece add +0.89 cm H₂O

Anesthetic gas scavenging

Type	Hospital system required	Machine connection
Active low flow:	High vacuum 36 L/min (300 mmHg) @ 12 in Hg	DISS evac
Active low flow:	Adjustable Venturi with flowmeter > 30 L/min	12.7 mm/ 0.5 in hose barb
Active high flow:	Low vacuum 40 - 130 L/min	30 mm/1.2 in BSI Male threaded
Active high flow:	Venturi 50 L/min	25 mm/0.98 in hose barb
Passive:	Passive or externally attached active system	30 mm/1.2 in M ISO taper
Active:	Venturi/Ejector > 30 L/min	12 mm/0.47 in hose barb
Active:	Venturi/Ejector > 30 L/min	8 mm/0.31 in hose barb
Active adjustable flow:	> 30 L/min	



Datex-Ohmeda Division • Instrumentarium Corp.
P.O. Box 900, FIN-00031 Datex-Ohmeda, Finland
Tel. +358 10 39411 • Fax +358 9 146 3310

www.datex-ohmeda.com

